
AlaxoStent

BRUGSANVISNING
GEBRAUCHSANWEISUNG
INSTRUCTIONS FOR USE
MANUAL DE INSTRUCCIONES
KÄYTTÖOHJEET
MANUEL D'UTILISATION
ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
UPUTE ZA UPOTREBU
HASZNÁLATI ÚTMUTATÓ
ISTRUZIONI PER L'USO
GEBRUIKSAANWIJZING
BRUKERVEILEDNING
INSTRUÇÕES DE USO
BRUKSANVISNING
KULLANMA TALIMATLARI



DA
DE
EN
ES
FI
FR
GR
HR
HU
IT
NL
NO
PT
SV
TR

REF

AP004 AlaxoStent C

service@alaxo.com
www.alaxo.com

AlaxoStent

- DA
- DE
- EN
- ES
- FI
- FR
- GR
- HR
- HU
- IT
- NL
- NO
- PT
- SV
- TR


DA	Komponenter af AlaxoStent
DE	Komponenten des AlaxoStent
EN	Components of AlaxoStent
ES	Componentes del AlaxoStent
FI	AlaxoStent-tuottoeen komponentit
FR	Composants de l'AlaxoStent
GR	Συστατικά μέρη της AlaxoStent
HR	Komponente proizvoda AlaxoStent
HU	Az AlaxoStent komponensei
IT	Componenti dello AlaxoStent
NL	Onderdelen van de AlaxoStent
NO	Komponentene til AlaxoStent
PT	Componentes de AlaxoStent
SV	AlaxoStent-systemets delar
TR	AlaxoStent Bileşenleri

DA	Sådan sættes de to dele af klipselukningen sammen
DE	Ineinanderschieben des Clipverschlusses
EN	Fitting into each other of the clip lock
ES	Acople del cierre de clip
FI	Pihtilukon osien sovittaminen toisiinsa
FR	Télescopage des deux parties du clip de fermeture
GR	Σύνδεση των δύο μερών με το κλιπ+
HR	Umetanje dijelova kopče za spajanje jedan u drugi
HU	a klipzár összekapcsolása
IT	Collegare la chiusura a clip
NL	In elkaar schuiven van de clipsluiting
NO	Delene av klemlåsen settes i hverandre
PT	Encaixando-se mutuamente no clipe de bloqueio
SV	Passa ihop klicklåsets ändar
TR	Klipsli kilidin birbirine takılması

DA	Nitinolgitret indsættes i indføringsrøret
DE	Einführen des Nitinolgeflechts in den Einführschlauch
EN	Insertion of the nitinol braid into the introduction tube
ES	Inserción de la malla de nitinol en el cateter de inserción
FI	Nitiinolitupoksen vieminen sisäänvientiputkeen
FR	Introduction du treillis de nitinol dans le tube d'introduction
GR	Εισαγωγή του πλέγματος νιτινόλης στο σωληνάκι εισαγωγής
HR	Uvođenje nitinolske mrežice u cijev vodilicu
HU	a nitinol szövét behúzása a bevezető tömlőbe
IT	Inserire la rete in nitinol nel tubo di inserzione
NL	Invoeren van het nitinolvlechtwerk in de invoerslang
NO	Stikk nitinolnettet inn i innføringslangen
PT	Inserção do entrançado nitinol no interior do tubo introdução
SV	För in nitinolnätet helt i införingslangen
TR	Nitinol ağ örgünün giriş tüpüne sokulması

DA
DE
EN
ES
FI
FR
GR
HR
HU
IT
NL
NO
PT
SV
TR

 4

DA
DE
EN
ES
FI
FR
GR
HR
HU
IT
NL
NO
PT
SV
TR

 5


DA	Illustration af placeringen af AlaxoStent i halsen; A=meatus, B=velum, C=tungen fundamentalt, D=epiglottis
DE	Schema für die Positionierung des AlaxoStent im Rachenraum; A=Nasengang, B=Gaumensegel, C=Zungengrund, D=Kehldeckel
EN	Scheme for positioning of the AlaxoStent in the throat; A=Nasal passage, B=Soft palate, C=Tongue base, D=Epiglottis
ES	Esquema de la colocación del AlaxoStent en la faringe; A=Meato nasal, B=Velo palatino, C=Base de la lengua, D=Epiglottis
FI	Scheme sijoittamiseksi AlaxoStent kurkussa; A=meatus, B=velum, C=kielen pohjimiltaan, D=kurkunkanteen
FR	Schéma pour le positionnement de l'AlaxoStent dans le pharynx; A=Conduit nasal, B=Voile du palais, C=Nase de la langue, D=Epiglote
GR	Σχήμα τοποθέτησης της Ενδοπρόθεσης AlaxoStent στο λαιμό; A=Ρινική δίοδος, B=Μαλακή υπερώα, C=Βάση γλώσσας, D= Επιγλωπίδα (Epiglottis)
HR	Prikaz postavljanja proizvoda AlaxoStent u grlo; A=nosnog hodnika, B=mekano nepce, C=jezik razlog, D=epiglotis
HU	Az AlaxoStent garatban történő elhelyezésének sémája; A=Orrjárat, B=Szájpadlás, C=Nyelvgyök, D= Gégefő (Epiglottis)
IT	Schema per il posizionamento di AlaxoStent nello spazio faringeo; A=Meato nasal, B=Velo palatino, C=Base della lingua, D=Epiglottide
NL	Schema voor de positionering van de AlaxoStent in de keelholte; A=Neusgang, B=Zachte gehemelte, C=Radix linguae, D=Strotklepje (Epiglottis)
NO	Plan for plassering av AlaxoStent i halsen; A=meatus, B=gane seil, C=tungen fundamentalt, D=epiglottis
PT	Esquema para o posicionamento do AlaxoStent na garganta; A=Passagem nasal, B=Palato mole, C=Base de língua, D=Epiglote
SV	Schematisk bild över placeringen av AlaxoStent i svalget; A=meatus, B=gommen segel, C=tungan bas, D=struplocket
TR	AlaxoStent'in boğazda konumlandırılma şeması; A=Burun pasajı, B=Yumuşak damak, C=Dil kökü, D=Gırtlak kapağı (Epiglottis)


DA	Fastsættelse af AlaxoStent med fikseringstape
DE	Fixierung des AlaxoStent mit dem Fixierpflaster
EN	Fixation of the AlaxoStent with the fixation tape
ES	Fijación de la AlaxoStent con la cinta fijadora
FI	Kiinnitys AlaxoStent kanssa kiinnitysteippi
FR	Fixation de LAlaxoStent avec le sparadrap fixatif
GR	στερέωση AlaxoStent με το ταινία στερέωσης
HR	Učvršćivanje AlaxoStent sa traka za pričvršćivanje
HU	Rögzítése AlaxoStent a rögzítő tapasz
IT	Fissaggio della AlaxoStent con il cerotto di fissaggio
NL	Bevestiging van de AlaxoStent met fixeerpleister
NO	Fixing AlaxoStent med festeteip
PT	Fixação do AlaxoStent com a fita de fixação
SV	Fastställande av AlaxoStent med fixeringstejp
TR	AlaxoStent'in sabitleme bandıyla sabitlemesi

DA
DE
EN
ES
FI
FR
GR
HR
HU
IT
NL
NO
PT
SV
TR

 6

 7

DA
DE
EN
ES
FI
FR
GR
HR
HU
IT
NL
NO
PT
SV
TR


DA	Nitinolgitret skubbes over rengøringsrøret
DE	Aufbringen des Nitinolgeflechts auf den Reinigungsschlauch
EN	Pushing of the nitinol braid over the cleaning tube
ES	Aplicación de la malla de nitinol en la tubo de limpieza
FI	Nitinolipunoksen työntäminen puhdistusputken päältä
FR	Mettre le treillis de nitinol sur le tuyau de nettoyage
GR	Σπρώχνουμε το πλέγμα πάνω στο καθαριστικό σωληνάκι
HR	Guranje nitinolske mrežice preko cijevi za čišćenje
HU	a nitinol szövet felhelyezése a tisztító tömlőre
IT	Inserimento della rete in nitinol sul tubo per la pulizia
NL	Opbrengen van het nitinolweefsel op de reinigings slang
NO	Nitinolnettet skyves over renseslangen
PT	Empurrando o entrançado nitinol sobre o tubo de limpeza
SV	Nitinolnätet träs över rengöringslangen
TR	Nitinol ağ örgünün temizleme tüpü üzerine itilmesi


DA	Nitinolgitret skubbes af rengøringsrøret
DE	Herunterschieben des Nitinolgeflechts vom Reinigungsschlauch
EN	Pushing off of the nitinol braid from the cleaning tube
ES	Deslizamiento de la malla de nitinol sobre el tubo de limpieza
FI	Nitinolipunoksen työntäminen puhdistusputkelta
FR	Retrait du treillis de nitinol du tuyau de nettoyage
GR	Αφαίρεση του πλέγματος νιτινόλης από το καθαριστικό σωληνάκι
HR	Guranje nitinolske mrežice s cijevi za čišćenje
HU	A nitinol szövet letolása tisztító tömlőről
IT	Spinta della rete di nitinol dal tubo per la pulizia
NL	Het van de reinigings slang afschuiven van het nitinolvechtwerk
NO	Nitinolnettet skyves av fra renseslangen
PT	Retirando o entrançado de nitinol do tubo de limpeza
SV	Nitinolnätet tas bort från rengöringslangen
TR	Nitinol ağ örgünün temizleme tüpünden itilmesi

DA
DE
EN
ES
FI
FR
GR
HR
HU
IT
NL
NO
PT
SV
TR

DA
DE
EN
ES
FI
FR
GR
HR
HU
IT
NL
NO
PT
SV
TR

DA	Symboler
DE	Symbole
EN	Symbols
ES	Simbolos
FI	Symbolit
FR	Symboles
GR	σύμβολα
HR	Simboli
HU	Szimbólumok
IT	Simboli
NL	Symbolen
NO	Symboler
PT	Simbolos
SV	Symboler
TR	Semboller

DA	Se brugsanvisningen
DE	Gebrauchsanweisung beachten
EN	Consult instructions for use
ES	Atender el manual de instrucciones
FI	Katso käyttöohjeet
FR	Respecter le manuel d'utilisation
GR	Συμβουλευτείτε τις οδηγίες χρήσης
HR	Pročitajte upute prije upotrebe
HU	A használati útmutató figyelembevételre
IT	Leggere attentamente le istruzioni per l'uso
NL	Neem de gebruiksaanwijzing in acht
NO	Se brukerveiledningen
PT	Consultar as instruções de uso
SV	Se bruksanvisningen
TR	Kullanma talimatlarına bakın


DA	Batchkode
DE	Chargennummer
EN	Batch code
ES	Número de lote
FI	Eräkoodi
FR	Numéro de lot
GR	Κωδικός παρτίδας
HR	Šifra serije
HU	Tételszám
IT	Lotto nr.
NL	Lotnummer
NO	Batchkode
PT	Código do lote
SV	Satskod
TR	Parti kodu

DA	Ordrenummer
DE	Bestellnummer
EN	Order number
ES	Número de pedido
FI	Tilausnumero
FR	Numéro de commande
GR	Αριθμός Παραγγελίας
HR	Broj za narudžbu
HU	Megrendelési szám
IT	Numero d'ordine
NL	Bestelnummer
NO	Ordrenummer
PT	Número da encomenda
SV	Beställningsnummer
TR	Sipariş numarası

DA
DE
EN
ES
FI
FR
GR
HR
HU
IT
NL
NO
PT
SV
TR

DA	lagring ved 4° til 50°C
DE	Lagerung bei 4° bis 50°C
EN	Temp. limitation: store between 4° and 50°C
ES	Almacenamiento a temperatura de 4° a 50° C
FI	Lämpötilarajoitus: Säilytä 4 ja 50°C välillä
FR	Conservé à une température entre 4° et 50°C
GR	να φυλάσσεται μεταξύ 4° και 50°C
HR	Temperaturno ograničenje: čuvati između 4° i 50°C
HU	Tárolás 4° és 50°C hőmérséklet között
IT	Conservare a una temperatura da 4° a 50°C
NL	Bewaren bij een temperatuur van 4° tot 50°
NO	Temperaturbegrensning: Oppbevares mellom 4°C og 50°C
PT	Limite de temperatura: armazenar entre 4° e 50°C
SV	Temperaturbegränsning: Förvara vid mellan 4°C och 50°C
TR	Sıcaklık sınırı: 4° ve 50°C arasında saklayın

DA

DE

EN

ES

FI

FR

GR

HR

HU

IT

NL

NO

PT

SV

TR

DA	Produktionsdato: Måned/år
DE	Herstellungsdatum: Monat/Jahr
EN	Manufacturing date: Month/Year
ES	Fecha de fabricación: mes/año
FI	Valmistuspäivämäärä: Kuukausi/vuosi
FR	Date de fabrication: mois/année
GR	ημερομηνία κατασκευής: Μήνας/Έτος
HR	Datum proizvodnje: mjesec/godina
HU	Gyártás időpontja: év/hónap
IT	Data di produzione: mese/anno
NL	Fabricagedatum: maand/jaar
NO	Produksjonsdato: Måned/år
PT	Data de produção: Mês/Ano
SV	Tillverkningsdatum: månad/år
TR	Üretim tarihi: Ay/Yıl

DA	Anvendes inden: Måned/år
DE	Verwendbar bis: Monat/Jahr
EN	Use by: Month/Year
ES	CAD: mes/año
FI	Käytettävä ennen: Kuukausi/vuosi
FR	À utiliser avant le: mois/année
GR	Χρήση έως: Μήνας/Έτος
HR	Upotrijebiti do: mjesec/godina
HU	Felhasználható: év/hónap
IT	Scadenza: mese/anno
NL	Te gebruiken tot: maand/jaar
NO	Bruk innen: Måned/år
PT	Uso por: Mês/Ano
SV	Använd före: månad/år
TR	Son kullanim: Ay/Yıl

DA	CE-mærke
DE	CE Kennzeichnung
EN	CE Mark
ES	Marca CE
FI	CE-merkintä
FR	Marquage CE
GR	Σήμα CE
HR	Oznaka CE
HU	CE jelzés
IT	Marcatura CE
NL	CE-markering
NO	CE-merke
PT	Marca CE
SV	CE-märkning
TR	CE işareti

DA	Producent
DE	Hersteller
EN	Manufacturer
ES	Fabricante
FI	Valmistaja
FR	Fabricant
GR	Κατασκευαστής
HR	Proizvođač
HU	Gyártó
IT	Produttore
NL	Fabrikant
NO	Produsent
PT	Fabricante
SV	Tillverkare
TR	Üretici

DA	Opbevares utilgængeligt for børn
DE	Für Kinder unzugänglich aufbewahren
EN	Store inaccessible for children
ES	Mantener fuera del alcance de los niños
FI	Säilytä lasten ulottumattomissa
FR	Tenir hors de la portée des enfants
GR	Να φυλάσσεται μακριά από παιδιά
HR	Čuvati izvan dohvata djece
HU	Gyermekek elől elzárva tárolandó
IT	Tenere fuori dalla portata dei bambini
NL	Buiten bereik van kinderen bewaren
NO	Oppbevares utilgjengelig for barn
PT	Armazenamento inacessível para as crianças
SV	Förvara oåtkomligt för barn
TR	Çocukların ulaşamayacağı şekilde saklayın

DA

DE

EN

ES

FI

FR

GR

HR

HU

IT

NL

NO

PT

SV

TR

1 Anvisninger

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INSTRUCTIONS FOR USE

1 Instructions

The AlaxoStent is a CE marked prescription medical device. It must only be used according to the briefing of the attending doctor.

1.1 CE Mark

The product bears the CE mark. The conformity evaluation assures that the product complies with the requirements of the EU directive for medical devices 93/42/EWG. According to the classification rules in Annex IX of the EU directive the AlaxoStent is a Class I product.

1.2 General Instructions

- The instructions for use is a component of the product. It has to be accessible to the user near to the product. Close adherence to these instructions is precondition for intended use and correct handling of the product.
- In case of damage resulting from use of third-party components or consumables no warranty claims apply.
- Alaxo is only liable for the product in respect of safety, reliability and functionality if it is used in compliance with the instructions for use.
- The instructions for use corresponds to the design of the product and the applicable safety norms as of going to print. All rights are reserved for the specified parts and components, processes, names and products.
- Full or partial reproduction subject to prior written approval by Alaxo.

1.3 General safety instructions

- The product has been developed and designed by Alaxo in such a manner that hazards during intended use are eliminated as far as possible. Nevertheless, we strongly recommend to strictly adhere to the following safety precautions in order to further reduce any residual risks to a minimum.
- It is not allowed to modify or reconstruct the product. In that case the CE mark lapses. Use of altered products is prohibited. In addition, Alaxo cannot take any warranty or liability for modified or reconstructed products. In order to assure safe application of the product users are liable for adherence to the instructions.

DA
DE
EN
ES
FI
FR
GR
HR
HU
IT
NL
NO
PT
SV
TR

- The original packaging should be stored for any potential return shipment. Only the original packaging ensures optimal protection of the product during transport. Should return shipment become necessary during the warranty period Alaxo does not take any liability for damages resulting from insufficient packaging.
- The product must not be accessible to children.
- The user has to inspect the product for functional safety and good condition each time prior to application.
- The user has to be familiar with handling of the product.
- The cleaning instructions have to be observed and obeyed to. Other cleaning procedures must not be applied. The cleaning spray (Prontosan® Wound Spray) must only be used for disinfectant cleaning of the AlaxoStent and the accessories. Do not swallow!
- The AlaxoStent must be used only by patients of age of 18 years and higher.
- The AlaxoStent is not suitable for life-supporting application in emergency care.

2.2 Contraindications

- Inappropriate structure of the nasal passages, especially if leading to problems during insertion of the introduction tube into a nasal passage
- Strong deviations of the nasal septum
- Dry rhinitis
- Disposition for nose hemorrhage
- Empty Nose Syndrome
- Shall not be used by persons with reduced mental ability, with reduced motor skills, with psychic restrictions, mental disorders or spastic syndromes

1.4 Symbols and indications

- One Prontosan® Wound Spray bottle may be used for 12 months after first use.
- The AlaxoStent may be used for max. 18 months after the first use.
- In the instructions for use, on the packaging and on the product the figures and symbols in  are used.

1.5 Precautionary measures

The attending doctor has to explain handling of the AlaxoStent to the patient. Especially, the doctor has to explain correct application and positioning of the AlaxoStent in the throat and its removal. The attending doctor has to advise on potential risks and complications connected to the AlaxoStent therapy and has to evaluate if the AlaxoStent therapy is applicable to the specific patient.

2 Product information

2.1 Intended use

- The AlaxoStent is a medical device and self-applied by the patient before going to sleep. The AlaxoStent is introduced via a nasal passage into the throat and positioned there in order to keep open the airway during sleep by mechanical splinting.
- The AlaxoStent is a prescription medical device. The attending doctor has to be sufficiently qualified in treatment of sleep apnea.
- Only one AlaxoStent is introduced through one nasal passage, never two or more at the same time.
- The AlaxoStent is clinically tested and CE marked for treatment of patients with obstructive sleep apnea or snoring.
- The attending doctor has to determine the disease state of the patient and if the AlaxoStent may be a useful therapy for this patient. The physician further has to examine if the nasal passages of the patient allow introduction of the AlaxoStent. The physician shall inform the patient which of the two nasal passages (left or right) is preferred for introduction of the AlaxoStent.
- In case of known lung diseases a pneumologic medical examination is required to evaluate if the AlaxoStent therapy can be applied to the specific patient.
- The AlaxoStent must be used only by one single patient.
- The AlaxoStent is suitable for repeated use up to 18 hours per application.

2.3 Warnings

- Persons with allergic reactions to nickel titanium (nitinol) may suffer an allergic response to this stent.
- Persons with allergic reactions to polyhexanide or undecylenamidopropylbetaine may suffer an allergic response to the cleaning spray.
- Do not use the AlaxoStent if this leads to pain! Do not terminate therapy without consulting your attending doctor.

2.4 Unintended use

- The only permitted use is splinting of the airway by introduction through a nasal passage into the throat. Other uses are prohibited.
- The AlaxoStent must not be introduced into orifices of the body other than the nose (e.g. ears or mouth).
- The AlaxoStent braid must not be pulled over any body parts or limbs.
- The product may only be used for the approved medical purposes. It must not be used for any other purposes, especially non-medical purposes.

2.5 Product description

The AlaxoStent comprises the following components (see  1):

- 1) Self-expanding nitinol braid [1] with clip lock [1a] (nitinol is a shape-memory metal, a nickel-titanium alloy widely used for medical devices)
- 2) Pusher [2] with clip lock (the end with the two holes) [2a]
- 3) Introduction tube [3]
- 4) Cleaning tube [4]
- 5) Fixation tape
- 6) 0.1% Polyhexanide / 0.1% betaine cleaning spray (Prontosan® Wound Spray)

3 Scope of Delivery


- | | |
|---|------------------------|
| 1 Nitinol braid | 1 Pusher |
| 1 Introduction tube | 1 Cleaning tube |
| 1 Fixation tape | 1 Instructions for use |
| 1 Bottle cleaning spray (Prontosan® Wound Spray 75ml) | |

DA
DE
EN
ES
FI
FR
GR
HR
HU
IT
NL
NO
PT
SV
TR



4 Use

The AlaxoStent is a medical device and has been developed and clinically tested for self-application by the patient. The nitinol braid consists of thin and sensitive metal wires so that the braid needs to be treated carefully and gently. The following instructions for handling and cleaning have to be obeyed in order to maintain the nitinol braid in good order over the allowed period of use of max. 18 months.

4.1 Preparation of the AlaxoStent for introduction

1. The nitinol braid is connected to the pusher [2] by fitting into each other the two parts [1a] and [2a] of the clip lock (see  2).
2. The other end of the pusher [2b] is inserted into black end [3a] of the introduction tube [3] and the pusher is moved through the introduction tube. Then the nitinol braid is carefully and completely pulled into the introduction tube using the pusher whereby it self-compresses. The clip lock shall be located in the introduction tube and the end of the stent shortly before the black end. (see  3)

4.2 Introduction of the AlaxoStent into the throat



1. The introduction tube is carefully introduced with the black end ahead into the right or left nostril and advanced horizontally through the nasal passage according to  4 down to the throat (tilt the head back!). Exact positioning of the nitinol braid is explained to the patient by the prescribing doctor. In each case the AlaxoStent stays away from the epiglottis.
2. The smooth material and the bent form of the introduction tube serves for good steerability at the end of the nasal passage down to the throat and is optimally adapted to the human anatomy. Nevertheless, differences in simplicity of the introduction process may arise e.g. from anatomic variability from patient to patient.
3. In case of dry mucosal tissue in the nose pushing forward of the introduction tube may lead to increased friction. To wetten the nasal mucosa, after consultation of the attending physician, a nasal spray containing sea salt, dexpanthenol and/or hyaluronic acid may be applied horizontally into the nasal passage. The introduction tube smoothly moves forward along the wet mucosal tissue with low friction.
4. After reaching the right position the pusher is held with the one hand and the introduction tube is slowly retracted over the pusher with the other hand. The nitinol braid self-expands first in the throat and then in the nasal passage and attaches itself towards the walls of the throat.
5. Once the introduction tube is fully removed the clip lock is pulled apart and thereby the pusher is separated from the nitinol braid. The AlaxoStent thus is in its position of use (see  4). Cleaning of the introduction tube is described in Section 4.4.
6. To secure the position of the AlaxoStent in the throat the part of the nitinol braid protruding from the nose is affixed with the fixation tape at the upper lip, the nose or the cheek (see  5). It has to be made sure that the fixation tape securely adheres to the nitinol braid.
7. Application of the AlaxoStent is completed and the patient can go to sleep.

8. Due to the opening force of the nitinol braid the throat is kept open and a collapse of the airway is prevented. Differences in the efficacy of splinting of the airway may be observed from patient to patient.

4.3 Removal of the AlaxoStent

After sleep the fixation tape is cautiously removed and the nitinol braid retracted from the nose and throat.

4.4 Cleaning of the AlaxoStent

1. After retraction from the nose the **introduction tube** is rinsed in lukewarm (**not hot !!!**) flowing drinking water until all potentially adhering mucus is removed.
2. The outer surface of the introduction tube is sprayed with the cleaning spray (Prontosan® Wound Spray) which is also sprayed into the end of the introduction tube. The cleaning spray has to act for one minute and then the introduction tube is let dry.
3. The **nitinol braid** is rinsed in lukewarm (**not hot !!!**) flowing drinking water immediately after retraction from the nose until all potentially adhering mucus is mostly removed.
4. The nitinol braid is pushed over the cleaning tube [4] and potentially adhering contaminations are carefully removed in lukewarm flowing drinking water with the fingers. For bringing the nitinol braid onto the cleaning tube its open end is cautiously pushed over one end of the tube (**do not pull!!!!**). Subsequently push the whole nitinol braid over the tube as shown in  6. Removal of contaminations from the nitinol braid with the fingers has to be done **only in longitudinal direction, never in cross direction (i.e. not around the cleaning tube)** as otherwise the sensitive nitinol braid may be damaged.
5. Once the nitinol braid is completely cleaned it is cautiously pushed off the cleaning tube with the thumb and one finger without exerting any pressure from the lower end to the upper end (see  7). Please **do not pull** the closed end nor the clip lock of the nitinol braid as it will contract due to its design and cannot be pulled off the cleaning tube!
6. Subsequently the nitinol braid is sprayed with the cleaning spray (Prontosan® Wound Spray) which has to act for one minute. Then the nitinol braid is rinsed once more in lukewarm flowing drinking water, again sprayed with the cleaning spray and let dry.
7. The **cleaning tube** and the **pusher** are rinsed, too, in lukewarm flowing drinking water, sprayed with the cleaning spray (Prontosan® Wound Spray) which has to act for one minute and then let dry.

The three-step cleaning procedure for the AlaxoStent serves for disinfection. By extensive rinsing the contaminations are efficiently removed. The 0.1% polyhexanide / 0.1% betaine cleaning spray (Prontosan® Wound Spray) additionally reaches all interspaces due to its modified surface tension. The ingredients possess excellent mucosa compatibility and a broad pathogen spectrum.

DA
DE
EN
ES
FI
FR
GR
HR
HU
IT
NL
NO
PT
SV
TR

4.5 Storage of the AlaxoStent

1. The dried nitinol braid shall not be stored in the introduction tube to prevent it from unnecessary pressure load.
2. The AlaxoStent shall not be stored in direct sunlight and at temperatures below 4°C or above 40°C.

4.6 Period of use

1. The AlaxoStent after first use may be used for max. 18 months for nightly application.
2. After this period an increased risk of breakage of a nitinol wire due to material wear exists. Therefore, usage is no longer allowed.
3. The nitinol braid must be inspected for any damage prior to each use. If breakage of a wire is observed the nitinol braid must not be used anymore. Inspection is easiest performed on the cleaning tube during the cleaning procedure.
4. A cleaning spray bottle (Prontosan® Wound Spray) may be used for 12 months after first use. Expiry date see print on label.
5. After expiry of these allowed periods of use a new AlaxoStent has to be employed.

5 Disposal

The AlaxoStent can be disposed of in domestic waste.

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Alaxo GmbH