
# AlaxoStent

BRUGSANVISNING  
GEBRAUCHSANWEISUNG  
INSTRUCTIONS FOR USE  
MANUAL DE INSTRUCCIONES  
KÄYTTÖOHJEET  
MANUEL D'UTILISATION  
ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ  
UPUTE ZA UPOTREBU  
HASZNÁLATI ÚTMUTATÓ  
ISTRUZIONI PER L'USO  
GEBRUIKSAANWIJZING  
BRUKERVEILEDNING  
INSTRUÇÕES DE USO  
BRUKSANVISNING  
KULLANMA TALIMATLARI

DA  
DE  
EN  
ES  
FI  
FR  
GR  
HR  
HU  
IT  
NL  
NO  
PT  
SV  
TR

**REF**

AP004 AlaxoStent C

service@alaxo.com  
www.alaxo.com

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INSTRUCTIONS FOR USE

## 1 Instructions

The AlaxoStent is a CE marked prescription medical device. It must only be used according to the briefing of the attending doctor.

### 1.1 CE Mark

The product bears the CE mark. The conformity assessment ensures that the product complies with the requirements of the Regulation (EU) 2017/745 (MDR). According to the classification rules in Annex VIII of the MDR the AlaxoStent is a Class I product.

### 1.2 General Instructions

The instructions for use is a component of the product. It has to be accessible to the user near to the product. Close adherence to these instructions is precondition for intended use and correct handling of the product. In case of damage resulting from use of third-party components or consumables no warranty claims apply. Alaxo is only liable for the product in respect of safety, reliability and functionality if it is used in compliance with the instructions for use.

The instructions for use corresponds to the design of the product and the applicable safety norms as of going to print. All rights are reserved for the specified parts and components, processes, names and products.

- Full or partial reproduction subject to prior written approval by Alaxo.

### 1.3 General safety instructions

- The product has been developed and designed by Alaxo in such a manner that hazards during intended use are eliminated as far as possible. Nevertheless, we strongly recommend to strictly adhere to the following safety precautions in order to further reduce any residual risks to a minimum.
- It is not allowed to modify or reconstruct the product. In that case the CE mark lapses. Use of altered products is prohibited. In addition, Alaxo cannot take any warranty or liability for modified or reconstructed products. In order to assure safe application of the product users are liable for adherence to the instructions.
- The original packaging should be stored for any potential return shipment. Only the original packaging ensures optimal protection of the product during transport. Should return shipment become necessary during the warranty period Alaxo does not take any liability for damages resulting from insufficient packaging.
- The product must not be accessible to children.
- The user has to inspect the product for functional safety and good condition each time prior to application.
- The user has to be familiar with handling of the product.
- The cleaning instructions have to be observed and obeyed to. Other cleaning procedures must not be applied.

DA
DE
EN
ES
FI
FR
GR
HR
HU
IT
NL
NO
PT
SV
TR

## 1.4 Symbols and indications

- The AlaxoStent may be used for max. 18 months after the first use.
- In the instructions for use, on the packaging and on the product the figures and symbols in  are used.

## 1.5 Precautionary measures

The attending doctor has to explain handling of the AlaxoStent to the patient. Especially, the doctor has to explain correct application and positioning of the AlaxoStent in the throat and its removal. The attending doctor has to advise on potential risks and complications connected to the AlaxoStent therapy and has to evaluate if the AlaxoStent therapy is applicable to the specific patient.

## 2 Product information

### 2.1 Intended use

- The AlaxoStent is a medical device and self-applied by the patient before going to sleep. The AlaxoStent is introduced via a nasal passage into the throat and positioned there in order to keep open the airway during sleep by mechanical splinting.
- The AlaxoStent is a prescription medical device. The attending doctor has to be sufficiently qualified in treatment of sleep apnea.
- Only one AlaxoStent is introduced through one nasal passage, never two or more at the same time.
- The AlaxoStent is clinically tested and CE marked for treatment of patients with obstructive sleep apnea or snoring.
- The attending doctor has to determine the disease state of the patient and if the AlaxoStent may be a useful therapy for this patient. The physician further has to examine if the nasal passages of the patient allow introduction of the AlaxoStent. The physician shall inform the patient which of the two nasal passages (left or right) is preferred for introduction of the AlaxoStent.
- In case of known lung diseases a pneumologic medical examination is required to evaluate if the AlaxoStent therapy can be applied to the specific patient.
- The AlaxoStent must be used only by one single patient. The AlaxoStent is suitable for repeated use up to 18 hours per application.
- The AlaxoStent must be used only by patients of age of 18 years and higher.
- The AlaxoStent is not suitable for life-supporting application in emergency care.

## 2.2 Contraindications

- Inappropriate structure of the nasal passages, especially if leading to problems during insertion of the introduction tube into a nasal passage
- Strong deviations of the nasal septum
- Dry rhinitis
- Disposition for nose hemorrhage
- Empty Nose Syndrome
- Shall not be used by persons with reduced mental ability, with reduced motor skills, with psychic restrictions, mental disorders or spastic syndromes

## 2.3 Warnings

- Persons with allergic reactions to nickel titanium (nitinol) may suffer an allergic response to this stent.
- Do not use the AlaxoStent if this leads to pain! Do not terminate therapy without consulting your attending doctor.

## 2.4 Unintended use

- The only permitted use is splinting of the airway by introduction through a nasal passage into the throat. Other uses are prohibited.
- The AlaxoStent must not be introduced into orifices of the body other than the nose (e.g. ears or mouth).
- The AlaxoStent braid must not be pulled over any body parts or limbs.
- The product may only be used for the approved medical purposes. It must not be used for any other purposes, especially non-medical purposes.

## 2.5 Product description

The AlaxoStent comprises the following components (see  1):

- 1) Self-expanding nitinol braid [1] with clip lock [1a] (nitinol is a shape-memory metal, a nickel-titanium alloy widely used for medical devices)
- 2) Pusher [2] with clip lock (the end with the two holes) [2a]
- 3) Introduction tube [3]
- 4) Cleaning tube [4]
- 5) Fixation tape

## 3 Scope of Delivery

- |                     |                        |
|---------------------|------------------------|
| 1 Nitinol braid     | 1 Pusher               |
| 1 Introduction tube | 1 Cleaning tube        |
| 1 Fixation tape     | 1 Instructions for use |

DA

DE

EN

ES

FI

FR

GR

HR

HU

IT

NL

NO

PT


SV

TR

## 4 Use

The AlaxoStent is a medical device and has been developed and clinically tested for self-application by the patient. The nitinol braid consists of thin and sensitive metal wires so that the braid needs to be treated carefully and gently. The following instructions for handling and cleaning have to be obeyed to in order to maintain the nitinol braid in good order over the allowed period of use of max. 18 months.

### 4.1 Preparation of the AlaxoStent for introduction

1. The nitinol braid is connected to the pusher [2] by fitting into each other the two parts [1a] and [2a] of the clip lock (see  2).
2. The other end of the pusher [2b] is inserted into black end [3a] of the introduction tube [3] and the pusher is moved through the introduction tube. Then the nitinol braid is carefully and completely pulled into the introduction tube using the pusher whereby it self-compresses. The clip lock shall be located in the introduction tube and the end of the stent shortly before the black end. (see  3)

### 4.2 Introduction of the AlaxoStent into the throat

1. The introduction tube is carefully introduced with the black end ahead into the right or left nostril and advanced horizontally through the nasal passage according to  4 down to the throat (tilt the head back!). Exact positioning of the nitinol braid is explained to the patient by the prescribing doctor. In each case the AlaxoStent stays away from the epiglottis.
2. The smooth material and the bent form of the introduction tube serves for good steerability at the end of the nasal passage down to the throat and is optimally adapted to the human anatomy. Nevertheless, differences in simplicity of the introduction process may arise e.g. from anatomic variability from patient to patient.
3. In case of dry mucosal tissue in the nose pushing forward of the introduction tube may lead to increased friction. To wetten the nasal mucosa, after consultation of the attending physician, a nasal spray containing sea salt, dexpanthenol and/or hyaluronic acid may be applied horizontally into the nasal passage. The introduction tube smoothly moves forward along the wet mucosal tissue with low friction.
4. After reaching the right position the pusher is held with the one hand and the introduction tube is slowly retracted over the pusher with the other hand. The nitinol braid self-expands first in the throat and then in the nasal passage and attaches itself towards the walls of the throat.
5. Once the introduction tube is fully removed the clip lock is pulled apart and thereby the pusher is separated from the nitinol braid. The AlaxoStent thus is in its position of use (see  4). Cleaning of the introduction tube is described in Section 4.4.

6. To secure the position of the AlaxoStent in the throat the part of the nitinol braid protruding from the nose is affixed with the fixation tape at the upper lip, the nose or the cheek (see  5). It has to be made sure that the fixation tape securely adheres to the nitinol braid.
7. Application of the AlaxoStent is completed and the patient can go to sleep.
8. Due to the opening force of the nitinol braid the throat is kept open and a collapse of the airway is prevented. Differences in the efficacy of splinting of the airway may be observed from patient to patient.

### 4.3 Removal of the AlaxoStent

After sleep the fixation tape is cautiously removed and the nitinol braid retracted from the nose and throat.

### 4.4 Cleaning of the AlaxoStent

1. After retraction from the nose the **introduction tube** is rinsed in lukewarm (**not hot !!!**) flowing drinking water until all potentially adhering mucus is removed.
2. The **nitinol braid** is rinsed in lukewarm (**not hot !!!**) flowing drinking water immediately after retraction from the nose until all potentially adhering mucus is mostly removed.
3. The nitinol braid is pushed over the cleaning tube [4] and potentially adhering contaminations are carefully removed in lukewarm flowing drinking water with the fingers. For bringing the nitinol braid onto the cleaning tube its open end is cautiously pushed over one end of the tube (**do not pull!!!!**). Subsequently push the whole nitinol braid over the tube as shown in  6. Removal of contaminations from the nitinol braid with the fingers has to be done **only in longitudinal direction, never in cross direction (i.e. not around the cleaning tube)** as otherwise the sensitive nitinol braid may be damaged.
4. Once the nitinol braid is completely cleaned it is cautiously pushed off the cleaning tube with the thumb and one finger without exerting any pressure from the lower end to the upper end (see  7). Please **do not pull** the closed end nor the clip lock of the nitinol braid as it will contract due to its design and cannot be pulled off the cleaning tube and let dry.
5. The **cleaning tube** and the **pusher** are rinsed, too, in lukewarm flowing drinking water and then let dry.

### 4.5 Storage of the AlaxoStent

1. The dried nitinol braid shall not be stored in the introduction tube to prevent it from unnecessary pressure load.
2. The AlaxoStent shall not be stored in direct sunlight and at temperatures below 4°C or above 40°C.

DA
DE
EN
ES
FI
FR
GR
HR
HU
IT
NL
NO
PT
SV
TR

#### 4.5 Storage of the AlaxoStent

1. The dried nitinol braid shall not be stored in the introduction tube to prevent it from unnecessary pressure load.
2. The AlaxoStent shall not be stored in direct sunlight and at temperatures below 4°C or above 40°C.

#### 4.6 Period of use

1. The AlaxoStent after first use may be used for max. 18 months for nightly application.
2. After this period an increased risk of breakage of a nitinol wire due to material wear exists. Therefore, usage is no longer allowed.
3. The nitinol braid must be inspected for any damage prior to each use. If breakage of a wire is observed the nitinol braid must not be used anymore. Inspection is easiest performed on the cleaning tube during the cleaning procedure.
4. After expiry of these allowed periods of use a new AlaxoStent has to be employed.

#### 5 Disposal

The AlaxoStent can be disposed of in domestic waste.

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Alaxo GmbH